



Ellman<sup>®</sup> introduces

# pelleFirm<sup>®</sup>

RF BODY TREATMENT

The NEW PelleFirm system provides a non-invasive, pain free RF treatment for the body, resulting in deep heating and the impermanent reduction in the appearance of cellulite.

Powered by the Pellevé<sup>®</sup> S5 RF generator, PelleFirm is ideal for problem zones including the abdomen, arms, buttocks and thighs. PelleFirm is Ellman's latest innovation to help you provide patients with healthier looking bodies.

PelleFirm uses a combination of RF energy for deep tissue heating and mechanical massage for cellulite. After just 4-6 treatments, doctors report that patients show significant improvement in the appearance of cellulite.



ellman<sup>®</sup>  
A Cynosure<sup>®</sup> Company

# pelleFirm®

RF BODY TREATMENT



30 mm electrode tip

PelleFirm handpieces are ergonomically designed to fit easily in your hand. Smaller than competitive devices, they are lightweight and versatile, allowing customized treatments across all contours of the body.

The small gentle spheres of the massage head facilitate movement of fluids into the interstitial space where it can be naturally evacuated by the body's lymphatic system.

25 mm & 30 mm electrode tips deliver RF energy to gently and effectively heat the skin up to 45 Celsius.



25 mm electrode tip



Before



After 4 Treatments

Courtesy of Igor Jeremic, MD

*"I have seen PelleFirm results in as little as a few treatments. We are very satisfied users."*  
Josephine Hernandez, MD



Before



After 5 Treatments

Courtesy of Josephine Hernandez, MD

Code	Description
P3PFVP-WMHU*	PelleFirm Vari-Pak w/Massage Heads [2 HP (25 & 30 mm) & 2 MH (25 & 30 mm)]
P325PFHPU*	PelleFirm 25 mm Handpieces (3ea)
P330PFHPU*	PelleFirm 30 mm Handpieces (3ea)
P325PFMH	PelleFirm 25 mm Massage Head (1ea)
P330PFMH	PelleFirm 30 mm Massage Head (1ea)

**Indications for Use:**

PelleFirm System has been CE marked for distribution in the EU with the following indications:

The PelleFirm device is intended to provide heating for the purpose of elevating tissue temperature for the treatment of selected medical conditions such as relief of pain, muscle spasms, and increase in local circulation.

The PelleFirm device is intended to provide skin tightening to reduce skin laxity and wrinkles, as well as a temporary reduction in the appearance of cellulite.

PelleFirm has been FDA-cleared in the United States for the following indications:

The PelleFirm RF device is intended to provide heating for the purpose of elevating tissue temperature for selected medical conditions such as temporary relief of pain, muscle spasms, and increase in local circulation. The PelleFirm massage device is intended to provide a temporary reduction in the appearance of cellulite.

\*For orders outside the U.S. drop the "U" at the end of the code (e.g. use P325PFHP and not P325PFHPU)

ellman®  
A CYNOSURE® Company

400 Karin Lane, Hicksville, NY 11801

P 516.594.3333 | F 516.279.4483 | www.ellman.com

Cynosure, Ellman, Pellevé and PelleFirm are registered trademarks of Cynosure, Inc.

©2014 Cynosure, Inc.  
CC14041Y REV C